

**SPECIFICATIONS**

	C1	C5 <sup>2</sup>	C6 <sup>2</sup>	C-101	LC-101 Euro Leather	D30
<b>Back Width:</b>	51 cm (20")	71 cm (28")	71 cm (28")	72 cm (28.5")	58 cm (23")	74 cm (29")
<b>Seat-To-Floor:</b>	41 cm (16")	46 cm (18")	48 cm (19")	47 cm (18.5")	47 cm (18.5")	48 cm (19")
<b>Seat Depth:</b>	51 cm (20")	46 - 51 cm (18 - 20")	53 cm (21")	51 cm (20")	51 cm (20")	56 cm (22")
<b>Seat Width:</b>	46 cm (18")	56 cm (22")	56 cm (22")	51 cm (20")	51 cm (20") <sup>1</sup>	56 cm (22")
<b>Top of Back to Seat:</b>	66 cm (26")	66 - 71 cm (26 - 28")	66 - 71 cm (26-28")	70 cm (27.5")	70 cm (27.5")	66 cm (26")
<b>Weight Capacity:</b>	113 kg (250 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)	147 kg (325 lbs.)	147 kg (325 lbs.)	170 kg (375 lbs.)
<b>Colour Options:</b>	Arctic Blue Chocolate Oatmeal	Arctic Blue Chocolate Oatmeal	Arctic Blue Chocolate Oatmeal <b>Levante Vinyl:</b> Royal Blue	Cocoa Fern Sandal	<b>Leather:</b> Beige Brown Black Red	Ocean Burgundy Spruce Latte <b>Levante Vinyl:</b> Burgundy Royal Blue Pine Mountain

	560	660	LL-805	T-3	LC-440 Power Stand Assist
<b>Back Width:</b>	71 cm (28")	71 cm (28")	77.5 cm (30.5")	69 cm (27")	53.5 cm (21")
<b>Seat-To-Floor:</b>	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	49 cm (19")
<b>Seat Depth:</b>	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	46 cm (18")
<b>Seat Width:</b>	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	54.5 cm (21.5")
<b>Top of Back to Seat:</b>	63.5 cm (25")	63.5 cm (25")	68.5 cm (27")	68.5 cm (27")	69 cm (27")
<b>Weight Capacity:</b>	170 kg (375 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)
<b>Colour Options:</b>	Ocean Burgundy Spruce	Latte	Cocoa Fern Sandal	Black (Leather)	Cocoa Fern Sandal

**BENEFITS**

- Engineered furniture laminate/hardwood frame
- Patented, quiet & smooth lift system
- Large side pocket <sup>1</sup>
- Hand control can be easily moved to left-hand side
- Integrated emergency battery backup (batteries not included)
- External transformer for easy access to battery backup

Note: Measurements may vary +/- 2.5 cm (1"). Specifications subject to change without prior notice.

1. Large side pocket not available on LC-440.
2. C5 and C6 also available in Pride's custom fabric selections range. Please see "Custom Fabric Selections" brochure for further details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



**Mobility**  
Hire & Sales

3 Blackwall Rd, Woy Woy - Ph (02) 4342 53087 - Email: enquiries@mobilityhireandsales.com.au



**Pride® Lift Chairs**  
Experience the Luxury

# Pride Lift Chairs

*Experience the Luxury*

Pride® Lift Chairs are designed to be the ultimate in lift chair comfort, style and performance and are available in several beautiful fabrics and styles to fit any décor. The state-of-the-art lift system utilised in Pride Lift Chairs will recline you further back than any other lift chair. Plus, Pride Lift Chairs are equipped with an emergency battery backup system should the power fail in your home.

## LIFT CHAIR DESCRIPTIONS

**C1**  
3-position, petite to small chaise lounger

**C5**  
3-position, medium to large chaise lounger with fully removable upholstery components

**C6**  
Infinite-position, medium to large chaise lounger with fully removable upholstery components and dual motors which enable independent recline/footrest function

**C-101**  
3-position, medium to large lift chair

**LC-101 Euro Leather**  
3-position, medium to large lift chair

**D30**  
3-position, medium to large chaise lounger

**560**  
3-position, small to medium chaise lounger




**660**  
Infinite-position, small to medium chaise lounger with dual motors that enable independent recline/footrest function

**LL-805**  
2-position, medium to large wall hugger chaise lounger which requires only 10 cm (4") of space between wall and chair

**T-3**  
Infinite-position, medium to large lift chair featuring Trendelenburg position (elevates legs above head for greater relaxation) and zero gravity position

**LC-440 Power Stand Assist<sup>3</sup>**  
Power seat lift

## AVAILABLE FABRICS <sup>1</sup>



## LIFT CHAIR MODELS



**C1**  
Accented Pillow-back Design  
Shown in Arctic Blue



**C5**  
Sewn Pillow-back Design  
Shown in Arctic Blue



**C6**  
Sewn Pillow-back Design  
Shown in Chocolate



**C-101**  
Button-back Design  
Shown in Sandal



**LC-101 Euro Leather**  
Split-back Design  
Shown in Beige Leather



**D30**  
Split-back Design  
Shown in Burgundy



**560**  
Pillow-back Design  
Shown in Spruce



**660**  
Pillow-back Design  
Shown in Latte



**805**  
Sewn Split-back Design  
Shown in Cocoa



**T-3**  
Mongolian-back Design  
Shown in Black (Leather)



**LC-440 Power Stand Assist<sup>2</sup>**  
Shown in Sandal



**PLEASE NOTE:** 2-position is a TV recline.  
3-position is a full recline.  
Infinite-position is a full recline with independent leg rest function.

<sup>1</sup> The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

<sup>2</sup> Not all fabrics are available for every model. Please see Specifications for details.

<sup>3</sup> This chair does not recline.