

HOMECRAFT ATLANTIC WAVE DELUXE SHOWER COMMODE, 18" ATTENDANT PROPELLED, ADJUSTABLE HEIGHT



This deluxe shower commode chair features a lightweight, silver painted aluminium frame, which provides a rust free, high strength product. Brakes on all four wheels for security, the swing away, detachable footrests offer comfort and support whilst being transported and the armrests can be flipped back to facilitate side transfer. This chair can be used in a shower cubicle or wet room, used as a commode or over a standard toilet bowl. The front of the seat is cut-away to facilitate personal hygiene and a set insert is provided for comfort during transport. Supplied with a square commode pan as standard. The detachable back enables the shower commode chairs to be partially dismantled for more compact storage.

FEATURES

- Chair can be used as commode, shower chair, or elevated toilet seat
- Height adjustable from 958 to 1008mm (with castors).
- Brakes on all four wheels
- Swing away arms and detachable footrests make transfers easy
- Helpful bathroom aid for the elderly and people with disabilities, spare parts are available

SPECIFICATIONS

Overall Height:	958 – 1008mm (with castors)
Overall Length:	550mm (without footrests) 810mm (with footrests)
Overall Width:	555mm
Seat Height:	495 – 570mm
Seat Width:	457mm
Seat Depth:	425mm
Aperture Size:	190mm x 360mm

Backrest Width:	462mm
Backrest Height (from seat):	436mm
Front & Rear Wheel Diameter:	130mm
Under Seat Clearance:	430 – 505mm
Weight:	12.4kg
Maximum User Weight:	125kg
Warranty:	12 months
ARTG Id:	199705

PRODUCT CODES

NEW-091563451	Homecraft Atlantic Wave Deluxe Shower Commode, 18" Attendant Propelled, Adjustable Height, Blue
SP-751011	PU Seat without Insert, Blue
SP-751011-LID	PU Seat Insert, Blue
SP-751003	PU Arm Rest Pad with Screws, Blue
SP-751003B-L	Arm Rest Assembly, Left, Blue
SP-751003B-R	Arm Rest Assembly, Right, Blue
SP-751023	Footrest, Left, Grey

SP-751024	Footrest, Right, Grey
SP-751031	130cm (5") Front Castor
SP-761031	130cm (5") Rear Castor
SP-751001	Backrest with Screws, Blue
SP-751004	Backrest Locking Knob
SP-CM5102209	PU Seat Clip with Screws
SP-751014	Bucket without Lid, White
SP-CM5176007	Locking Bracket with Screw for Footrest

CALL 1300 473 422

sales@performhealth.com.au | www.performhealth.com.au