

icare IC333

The IC333 Bed is unmatched for function and durability, taking the hospital look out of your room. This four function bed, features dual lift motors totalling 12,000 newtons lifting power, unlimited under-bed clearance and eight heavy duty castors.



KEY FUNCTIONS:




Head Lift



Knee Break



Height Adjustable
22-66cm



Trendelenburg and
Reverse Trendelenburg

BED CODES & SIZES:



LONG SINGLE

Code: IC333LS
Mattress Surface Size: 2030 x 900
Base Size: 2150 x 960



KING SINGLE

Code: IC333KS
Mattress Surface Size: 2030 x 1070
Base Size: 2150 x 1120



LONG DOUBLE

Code: IC333D
Mattress Surface Size: 2030 x 1350
Base Size: 2150 x 1400
One hand control



QUEEN

Code: IC333Q
Mattress Surface Size: 2030 x 1500
Base Size: 2150 x 1580
One hand control



SPLIT KING

Code: IC333K
Mattress Surface Size: 2030 x 900 (per side)
Base Size: 2150 x 960 (per side)
Individual hand control per side
Requires two long single headboards


COLOUR CHOICES:




Onyx




Stone




Silver



Smoke



Vinyl Ebony



Vinyl Dove

Note: Always remember that TWO Long Single beds together make a King size bed.

KEY FEATURES:



200kg SWL on Long Single And King Single



Strong Breathable Mattress Platform



Available With Headboards



External Rechargeable Battery Backup



Lockable Heavy Duty Castors



300kg SWL on Double and Queen



End and Corner Mattress Retainer Bars



10 Year Warranty



Back Lit Wipeable Hand Piece



No Gap Between Mattresses on King Beds (2x long singles)



Unlimited Under Bed Clearance



Ultra Low Design Reducing Risk of Falls



Optional 3 sided surround for wheelchair transfers.

TECHNICAL SPECIFICATIONS:

- » Safe working load - 200kg on long single and king single sizes and 300kg for double and queen sizes. The safe working load includes any model of Icare mattresses.
- » Hand Piece - Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.
- » Fabric - Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties and waterproof.
- » Mattress Surface - Strong metal mesh under mattress allowing increased airflow
- » Actuators - 1x 6000 newton - 24 watt head adjust actuator. 1x 6000 newton - 24 watt knee adjust actuator. 2x 6000 newton - 24 watt height adjust actuators (2x 8000 newton on double and queen sizes)
- » No pull motors for added safety.
- » Wheels - 8x heavy duty swivel castors with one lockable castor on each corner.
- » Head lift angle 62° and knee break angle 42°
- » Max height - 66cm from floor to top of base
- » Min height - 22cm from floor to top of base

ACCESSORIES AVAILABLE:

- » High Side Rail
- » Low Side Rail
- » Bed Stick
- » Over Bed Pole
- » Full Length Side Rails
- » Headboard
- » Footboard



Available in one piece Double and Queen