



Days Steel Over Toilet Aid, with Seat and Splashguard



The Days Over Toilet Aid assists transfers on & off the toilet. This toilet surround frame is height adjustable, with a push button adjustment. Powder coated steel frame with non slip rubber feet. Comes complete with a seat and splashguard.

FEATURES

- Powder coated steel toilet surround
- Adjustable in height
- Moulded plastic armrests
- To suit various end users
- Non-marking rubber tips, slip resistant for added safety

SPECIFICATIONS

- Seat Width: 370mm
- Width Between the Arms: 450mm
- Height Adjustable: 440mm - 590mm
- Maximum User Weight: 125kg
- Complies with AS/NZS 3973:2009 & AS/NZS 4688.2:2000 & AS/NZS 4688.3:2000.

PRODUCT CODES

JAN-AC20

Days Steel Over Toilet Aid, with Seat and Splashguard

JAN-AC20SP-009

Splashguard for JAN-AC20 Days Over Toilet Aid

JAN-AC20SP-008

Toilet Seat for JAN-AC20 Days Over Toilet Aid