

PRIDE POWER RECLINER COLLECTION



MODEL NAME	C1	C5	C6	C-101	LC-101 Euro Leather	LC-107	D30	560	660	T-3	LC-358XL
POSITION	3-Position Lift Chair 163 cm & below (5'4" & below)	3-Position Chaise Lounger with fully removable upholstery components 163 cm-183 cm (5'4" - 6'0")	Infinite-Position Chaise Lounger with fully removable upholstery components 163cm-183 cm (5'4" - 6'0")	3-Position Chaise Lounger 163 cm-183 cm (5'4" - 6'0")	3-Position Chaise Lounger with European style rounded arm design 163 cm -183 cm (5'4" - 6'0")	Infinite-Position Chaise Lounger with extended length on footrest 163 cm -183 cm (5'4" - 6'0")	3-Position Chaise Lounger 163 cm-183 cm (5'4" - 6'0")	3-Position Chaise Lounger 163 cm-183cm (5'4" - 6'0")	Infinite-Position Chaise Lounger 163 cm-183cm (5'4" - 6'0")	Infinite-Position Chaise Lounger featuring Trendelenburg position (elevates legs above head for greater relaxation) & zero-gravity position 163 cm-183 cm (5'4" - 6'0")	3-Position Chaise Lounger 163 cm-183 cm (5'4" - 6'0")
WEIGHT CAPACITY	113 kg (250 lbs)	170 kg (375 lbs)	136 kg (300 lbs)	147 kg (325 lbs)	147 kg (325 lbs)	147 kg (325 lbs)	170 kg (375 lbs)	170 kg (375 lbs)	170 kg (375 lbs)	170 kg (375 lbs)	227 kg (500 lbs)
BACK WIDTH	51 cm (20")	71 cm (28")	71 cm (28")	72 cm (28.5")	58 cm (23")	74 cm (29")	74 cm (29")	71 cm (28")	71 cm (28")	69 cm (27")	91 cm (36")
SEAT-TO-FLOOR	41 cm (16")	46 cm (18")	48 cm (19")	47 cm (18.5")	47 cm (18.5")	51 cm (20")	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	51 cm (20")
SEAT DEPTH	51 cm (20")	46 cm-51 cm (18" - 20")*	53 cm (21")	51 cm (20")	51 cm (20")	53 cm (21")	56 cm (22")	48 cm (19")	51 cm (20")	51 cm (20")	56 cm (22")
SEAT WIDTH	46 cm (18")	56 cm (22")	56 cm (22")	51 cm (20")	51 cm (20")	51 cm (20")	56 cm (22")	48 cm (19")	48 cm (19")	51 cm (20")	66 cm (26")
TOP OF BACK TO SEAT	66 cm (26")	66 cm-71 cm (26" - 28")*	66 cm-71 cm (26" - 28")*	70 cm (27.5")	70 cm (27.5")	70 cm (27.5")	66 cm (26")	63 cm (25")	63 cm (25")	68.5 cm (27")	74 cm (29")
BACK STYLE	Blown Fibre Accented Pillow-Back	Blown Fibre Sewn Pillow-Back	Blown Fibre Sewn Pillow-Back	Blown Fibre Button-Back	Blown Fibre Split-Back	Blown Fibre Sewn Split-Back	Blown Fibre Sewn Split-Back	Blown Fibre Pillow-Back	Blown Fibre Pillow-Back	Blown Fibre Mongolian-Back	Blown Fibre Button-Back
FULLY PADDED CHAISE	Yes	Yes ²	Yes ²	Yes	Yes	Yes	Yes ¹	Yes ¹	Yes ¹	Yes	Yes
INTEGRATED BATTERY BACKUP³	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
RECLINE POSITIONS	3	3	Infinite	3	3	Infinite	3	3	Infinite	Infinite	3
STANDARD ARM COVERS	No	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes
STANDARD HEAD COVERS	No	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes
FABRIC SELECTION GUIDE	<p>FABRIC</p> <p>Arctic Blue Chocolate Faun</p> <p>LEVANTE (VINYL)</p> <p>Royal Blue</p>	<p>FABRIC</p> <p>Arctic Blue Chocolate Faun</p>	<p>FABRIC</p> <p>Arctic Blue Chocolate Faun</p> <p>LEVANTE (VINYL)</p> <p>Royal Blue</p>	<p>FABRIC</p> <p>Cocoa Fern Sandal</p>	<p>LEATHER</p> <p>Beige Brown Black Red</p>	<p>FABRIC</p> <p>Pecan Burgundy Ocean</p>	<p>FABRIC</p> <p>Burgundy Ocean</p> <p>LEVANTE (VINYL)</p> <p>Royal Blue</p>	<p>FABRIC</p> <p>Burgundy Ocean</p> <p>LEVANTE (VINYL)</p> <p>Royal Blue</p>	<p>FABRIC</p> <p>Latte</p>	<p>LEATHER</p> <p>Black</p> <p>FABRIC</p> <p>Latte Burgundy</p>	

3 - POSITION	The 3-position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing
TRUE INFINITE-POSITION	The true infinite-position recline is perfect for taking a nap. Infinite-position Pride® Power Lift Recliners enable the backrest and footrest to move independently from one another so that a variety of positions can be reached.
TRENDELENBURG	The Trendelenburg position, elevates one's feet higher than his or her heart for increased comfort and relaxation.
ZERO GRAVITY	The zero gravity position, which can be reached on infinite-position Pride® Power Lift Recliner models, positions one's thighs and lower legs higher than his or her torso



NOTE The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Measurements may vary +/- 2.5cm (1") Specifications subject to change without prior notice. ¹ Removable fully padded chaise. ² All upholstery compounds on C5 and C6 are fully removable. ³ Batteries are not included. *Frame Adjustment