



AJ0020

Overtilet Aid - Aluminium



Features & Benefits

- Provides a raised seat and support arms to assist toilet transfers
- Ideal for people with hip or knee injuries or those who require greater seat height
- Ergonomic angled arms assist transfers
- Height adjustable
- Removable toilet seat and splashguard
- Non-slip rubber tips
- Powder coated aluminium frame is lightweight and rust proof
- Tested and conforms to Australian Standards AS/NZS 3973: 2009

Also Available

- AZ1400 Lid, Clip-on

Specifications

Seat Height	Width Between Arms	Arm Height	Base Width	Base Depth	Max User Weight
470-620mm	445mm	230mm	550mm	610mm	150kg