



Bedside Commode



The Days Bedside Commode can be stacked up to 10 high. It has a clip on toilet seat, which clips on in seconds. There are no metal fittings on toilet seat. It is also height adjustable and easy to clean. It has a light blue deluxe padded overseat, as well as a removable round bowl with lid and handle. The lid eliminates spillages and reduces risks of cross infection. It also features a powder coated dark blue frame. This bedside commode looks very smart in a bedroom.

SPECIFICATIONS

• Overall Height:	820 – 970mm	• Weight:	8kg
• Overall Width:	610mm	• Max User Weight:	190kg
• Overall Length:	550mm	• Warranty:	12 months
• Seat Height:	460 – 610mm	• ARTG Id:	n/a
• Width between Arm:	504mm		

PRODUCT CODES

SHH-DAYCOM521A	Days Bedside Commode, 3/Ctn
MFI-AC90B	Spare Parts - Plastic Bowl with Plastic Lid and Metal Handle
NOV-CHA71237	Spare Parts - Padded Backrest
SHH-521ASA003	Spare Parts - Padded Seat
SHH-521APL002	Spare Parts - Toilet Seat, Black
SHH-521A-000003	Spare Parts - Extension Leg Clip/Pin, Grey
SHH-521A-000004	Spare Parts - Sleeve at End of Leg, Grey
SHH-521A-000005	Spare Parts - Rubber Tip, Grey

CALL 1300 473 422

sales@performhealth.com.au | www.performhealth.com.au