icare

The IC111 Bed is our low risk extra mobility option. These domestic looking bases offer head lift and knee break. Split queen and split king sizes are available with split mattresses with no gap and individual controls.

KEY FUNCTIONS:



Knee Break

BED CODES & SIZES:



HALF QUEEN

Code: IC111HQ Mattress Surface Size: 2030 x 750 Base Size: 2150 x 780 (as partner) Ideal for partner configurations and small rooms



LONG SINGLE

Code: IC111LS Mattress Surface Size: 2030 x 900 Base Size: 2150 x 960



KING SINGLE

Code: IC111KS Mattress Surface Size: 2030 x 1070 Base Size: 2150 x 1120



SPLIT QUEEN

Code: IC111Q

Mattress Surface Size: 2030 x 750 (per side) Base Size: 2150 x 790 (per side)

- Individual hand control per side }>
- }> Fits a one piece queen headboard




SPLIT KING Code: IC111K

Mattress Surface Size: 2030 x 900 (per side) Base Size: 2150 x 960 (per side)

- Individual hand control per side }>
- }> Requires two long single headboards

COLOUR CHOICES:



Onyx

Smoke

Silver



KEY FEATURES:



200kg Safe Working Load Per Side



Strong Breathable



End And Corner

Bars

Fixed Base

Height 37cm



Available with Headboards



5 Year Warranty



External Rechargeable Battery Backup



Back Lit

Wipeable Hand

Piece



icare

Lockable Heavy Duty Castors



No Gap Between Mattresses On Queen And King Beds

TECHNICAL SPECIFICATIONS:

Mattress Platform Mattress Retainer

- » Safe working load The safe working load is 200kg per side and includes any model of lcare mattresses.
- » Hand Piece Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.
- » Fabric Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties and waterproof.
- » Mattress Surface Strong metal mesh under mattress allowing increased airflow.
- » Actuators 1x 6000 newton 24 watt head adjust actuator. 1x 6000 newton 24 watt knee adjust actuator.
- » No pull motors for added safety.
- » Wheels large lockable heavy duty swivel castors on each corner.
- » Head Lift angle 62°
- » Knee break angle 42°
- » Fixed base height 37cm from floor to top of base.

ACCESSORIES AVAILABLE:

- » High Side Rail
- » Low Side Rail
- » Bed Stick
- » Over Bed Pole
- » Full Length Side Rails
- » Headboard
- » Footboard

