

ZERO TURN10



Shown in Candy Apple Red

International Seat Shown

181.4 kg
Weight Capacity

Up To
10 km/h
Maximum Speed

109.2 cm
Turning Radius

Up To
38.6 km
Range Per Charge*



Pride's iTurn Technology™

The **Zero Turn 10** (ZT10) boasts a sleek, contemporary design that can traverse the most beautiful, out-of-reach places that can also be used in your home, moving seamlessly through doors and down hallways. Combining powerful dual motors and two-wheel drive with revolutionary turning technology, the ZT10 will carry you to places unknown.

Specifications

Model Name	Zero Turn 10 (ZT10)
Weight Capacity	181.4 kg (400 lbs.)
Maximum Speed ¹	10 km/h (6.2 mph)
Ground Clearance ³	5 cm (2") at motor
Turning Radius ³	109.2 cm (43")
Overall Length ³	121.9 cm (48")
Overall Width ³	61.6 cm (24.25")
Seat-to-Ground Height ³	58.4 - 63.5 cm (23 - 25")
Seat-to-Deck Height ³	45.7 cm - 50.8 cm (18 - 20")
Front Tyres ³	22.8 cm (9") solid
Rear Tyres ³	27.3 cm (10.75") solid
Suspension	Front & rear (CTS Suspension)
Range Per Charge ^{1,2}	Up to 28.8 km at 181.4 kg Up to 38.6 km at 90.7 kg
Total Weight Without Batteries ⁵	82.46 kg
Weight of Heaviest Piece ³	34 kg (front section)
Standard Seating	Type: Foldable, compact Weight: 13.6 kg Material: Black vinyl Dimensions: 45.7 x 43.18 cm
Drive System	24V DC motor, rear-wheel drive, dual in-line
Dual Braking System	Regenerative & electromechanical
Battery Requirements ⁴	Size: (2) 40 AH Weight: 13.7 kg (each)
Battery Charger	3.5 amp, off-board
Warranty	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.
2. Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
4. Battery weight may vary based on battery manufacturer.
5. Includes standard seat weight.

DISCLAIMER: The intended use of this Pride Mobility Products® Scooter is to provide mobility to persons with full use of their upper and lower extremities and who are able to support themselves in an upright position. This product is primarily for indoor use and the outdoor capabilities are limited to hard, flat, and even surfaces such as concrete, black top, and other similar surfaces. While we understand that the scooter's outdoor capabilities are limited, when driving near uneven obstacles, follow the safety guidelines outlined within the safety section of this manual.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius, and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Refinement meets innovation!

- Pride's patent pending, intelligent-turning, iTurn Technology provides 4-wheel stability with 3-wheel manoeuvrability
- Pride's CTS Suspension
- Three-speed control lets you adjust your speed when moving from indoor to outdoor driving
- Easy access tiller-mounted charger port
- USB charger built into the tiller for convenient charging of smart phones and portable devices
- Integrated storage/cup holder in the tiller
- New, improved ambient voltmeter
- Delta tiller with wraparound handles
- Standard large basket
- Full LED light package includes LED headlights and rear lights with sequential turn signals and convenience light under tiller
- Charging port convenience light
- Pinchless hinge, molded plastic, black vinyl seat

Available Accessories

- Rear bag / walking stick holder
- Walker holder*
- Oxygen tank holder*
- Mobile phone holder
- Cup holder
- Saddle Bag

* Must have Pride rear bag to attach

Available Colour



Candy
Apple
Red



Quality
Endorsed
Company
100% 05120205
SAI Global



20-24 Apollo Drive, Hallam VIC 3803
ACN 088 609 661 ABN 66 088 609 661

AU_ZeroTurn10_24-08-20