

# icare IC100

The IC100 Bed is a non-function hospital grade bed, with a static mattress platform. This bed provides a strong base for a partner or parent to sleep alongside the client.



## BED CODES & SIZES:



### LONG SINGLE

Code: IC100LS

Mattress Surface Size: 2030 x 900

Base Size: 2150 x 960



### KING SINGLE

Code: IC100KS

Mattress Surface Size: 2030 x 1070

Base Size: 2150 x 1120


## COLOUR CHOICES:




Onyx



Stone



Silver



Smoke



Vinyl Ebony



Vinyl Dove

## KEY FEATURES:



200kg Safe Working Load



Fixed Base Height 37cm



Available with Headboards



Lockable Heavy Duty Castors



Strong Breathable Mattress Platform



End And Corner Mattress Retainer Bars



10 Year Warranty

## TECHNICAL SPECIFICATIONS:

- » Safe working load - The safe working load is 200kg includes any model of Icare mattresses.
- » Fabric - Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties and waterproof.
- » Mattress Surface - Strong metal mesh under mattress allowing increased airflow.
- » Wheels - large lockable heavy duty swivel castors on each corner.
- » Fixed base height - 37cm from floor to top of base.

## ACCESSORIES AVAILABLE:

- » High Side Rail
- » Low Side Rail
- » Bed Stick
- » Over Bed Pole
- » Full Length Side Rails
- » Headboard
- » Footboard