


P-16 Pro Lite



6 km/h

Max. Speed

15 km

Travel Distance

110 kg

Weight Capacity

16.5 kg

Total Unit Weight

Specifications:

Weight Capacity:	110 kgs
Max. Speed:	6 km/h
Travel Distance:	15 km
Total Unit Length ¹ :	95 cm
Total Unit Width:	61 cm
Total Unit Weight:	16.5 kgs
Seat Dimensions (D x W x H):	41 x 50 x 48 cm
Drive System:	2x120W
Battery Requirements:	24V 10AH Lithium
Battery Charger:	DC24V/2.0A


¹ Measurement taken without footrests. Footrests measure 39 cm in length.

Lexham Mobility introduces a new lightweight and durable magnesium alloy power wheelchair. The compact and powerful P-16 Pro Lite marks a significant shift in transport-friendly mobility, benefitting from an easy-fold design and a disassembled weight of just 14.8 kgs (without battery and legrests).

Wherever your journey takes you, the tough, portable, and stylish Pro Lite is sure to deliver a new level of convenience and functionality.



Easy-fold
Frame



Simple
Disassembly



Convenient
Carry Bag

Visit lexhammobility.com to learn more.

LEXHAM
mobility