

Jazzy®

1450

The Jazzy® 1450 offers heavy-duty construction and performance components including upgraded motors, a reinforced frame, 90-amp PG VR2 controller and speeds up to 7.2 km/h (4.5 mph). This bariatric power chair features a 272 kg (600 lbs.) weight capacity, front-wheel drive for superb indoor maneuverability in the tightest of spaces and large front drive wheel for extra absorption over obstacles. The Jazzy 1450 combines superb climbing capabilities with front anti-tips for transitioning over grade changes and handling various terrains easily.

*Jazzy® 1450
shown in Candy Apple Red
with high-back seat*



SPECIFICATIONS

Weight capacity	272 kg (600 lbs.)
Maximum speed	Up to 7.2 km/h (4.5 mph)
Ground clearance	3.81 cm (1.5")
Turning radius	67.3 cm (26.5")
Overall length	97.16 cm (38.25")
Overall width	74.3 cm (29.25")
Drive wheels	35.53 cm (14") solid
Front wheels	10.16 cm (4") solid anti-tips
Rear wheels	22.9 cm (9") castors
Suspension	Swing Arm
Drivetrain	Two-motor, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	90 amp, PG VR2
Battery charger	Off-board, 8 amp
Per-charge range	Up to 19.62 km (12.19 miles)
Battery requirements	(2) 12 volt, Group 24
Battery weight	23.6 kg (52 lbs.) each
Base weight	88.5 kg (195 lbs.)
Standard seat weight	23.4 kg (51.5 lbs.) (high-back)
Standard seat	61 cm x 61 cm (24" x 24")
Max. seat dimensions	81 cm x 71 cm (32" x 28")

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

FEATURES

- Black vinyl high-back seat with height-adjustable headrest
- Flat-free drive tyres
- Heavy-duty, base mounted, flip-up height adjustable foot platform
- Dual-layer powder coated frame for increased durability

OPTIONS & ACCESSORIES

- Cup holder
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Safety flag



COLOUR CHOICES



Candy
Apple Red



Viper
Blue

Pride[®]
Mobility Products
Australia Pty. Ltd.

20-24 Apollo Drive Hallam, Victoria 3803
ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au

