



# Etac Beauty body care products

Award-winning ergonomically designed body care

Beauty is a light weight range of ergonomically designed body care products, which are all carefully balanced to provide maximum power with minimal effort and strain on hands, arms and shoulders.

Each product has a unique shape that makes reaching even easier.

**etac**<sup>®</sup>

Creating Possibilities

# Body Washer and Back Washer

Beauty Body Washer or Beauty Back Washer are designed to make it easy to reach and clean the body. The wash cloths do not absorb water which make the products light and balanced even when wet.

The handles are oval and have a surface which enables a secure grip. The wash cloths are replaceable.



# Hair Washer

Beauty Hair Washer substitutes the fingers to apply shampoo and massage the scalp. The shaft and handle are of the same design as the body washer, making it easier to reach the head without strain on the hands, shoulders or arms.



Description	Length	Weight	Order Code
Back Washer	770 mm	203 g	80210080
Back Washer, Replacement cloth			80210104
Body Washer	380 mm	112 g	80210076
Body Washer, Replacement cloth			80210103
Hair Washer	360 mm	112 g	80210078

### Material

Polypropylene and polypropylene based TPE  
Wash cloth: Nylon with polythene inner lining

### Care

Beauty wash cloths (nylon fabric and inner lining) are machine washable at 60°C.



# Comb and Hairbrush

Beauty Comb and Beauty Hairbrush come in two lengths. The shorter version is recommended when the user can reach up to the ear. The longer version is suitable for users who can only reach to the shoulder. Both the comb and brush are curved to follow the shape of the head. This means that the handle can be held close to the body to provide high comfort and control when combing/brushing the hair.



Description	Length	Weight	Order Code
Hairbrush	290 mm	102 g	80210072
Hairbrush, Long	370 mm	112 g	80210074
Comb	300 mm	80 g	80210068
Comb, Long	370 mm	90 g	80210070


**Material**  
Polypropylene and polypropylene based TPE

# Multipurpose grip

Beauty Multipurpose grip rests securely in the whole hand and offers a variety of use. With a poor hand function it can be very difficult to use a nail file. By inserting a nail file into the narrow end, a pedicure can be done with ease. The wider end of the grip may be used for a tooth brush or shaver.



*The Multipurpose Grip has a soft smooth surface.*



Description	Length	Weight	Order Code
Multipurpose Grip includes 3 nail files	110 mm	45 g	80210082

**Material**  
Polypropylene and polypropylene based TPE  
Nail file: Emery board

# Beauty Kit



Beauty Kit	Order Code
1 Body Washer 1 Comb Long 1 Multipurpose Grip with 3 nail files	80210066



 **etac**<sup>®</sup>  
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”