

Arm / Leg Shearing Sleeve



Cotton/lycra fabric

Prevents shearing and/or minor cuts



Elastic strap holds sleeve in position

A cotton lycra sleeve designed to fit snugly over the forearm to prevent the skin on the arm and back of the hand from shearing and/or minor cuts.

The sleeve may extend beyond the elbow or be folded to make the sleeve shorter. The sleeve can also be used to prevent shearing on the shins and lower leg.

An elasticised strap fits over the palm to hold the sleeve in position over the back of the hand (or fits around the instep to hold the sleeve in position over the top of the foot). If extra protection is needed over a particular area, additional padding can be placed around the arm or leg with the sleeve used over the top to hold the padding in place.

- › Machine washable up to 80°C
- › One size fits all



Code	P192
Length	350 mm
Width	90 > 130 mm
Weight	50 g
Conforms to	ARTG 100917
Warranty	1 year