

Proven pressure relief utilising the unique Air Comfort Seating System.



Large rear castors



Brake system



Rear push handle

Improving your life everyday™

Sales, administration
and general enquiries
National 1300 783 783

Your Customer Service Team
customerservice@kcare.com.au

kcare.com.au



COMPACT MOBILE LIFT CHAIRS

Improving your
life everyday™

KCare
Healthcare Solutions

Compact Mobile Lift Chairs

Features

- Total surface air pressure system for long term pressure relief
- Removable seat and back cushions
- Electric lift and recline mechanism
- Front and rear castors
- Wall Saver mechanism
- Back-up battery capability
- Easy to clean, repair and maintain
- Available in CarrEx
- Custom colours available
- Available in single or twin motor



Removable seats and back cushions for easy cleaning



Total surface air pressure system



Electric lift and recline mechanism from full recline to assisted standing



Description	Overall Width	Overall Height	Overall Length	Seat Depth	Seat Width	Seat Height	Back Height	Max User Weight
Small								
AC59098 (Single Motor) AC59100 (Twin Motor)	750mm	1060mm	820mm	440mm	500mm	480mm	680mm	140kg
Medium								
AC59087 (Single Motor) AC59090 (Twin Motor)	750mm	1060mm	820mm	510mm	500mm	820mm	740mm	140kg
Large								
AC59099 (Single Motor) AC59101 (Twin Motor)	820mm	1060mm	890mm	510mm	600mm	520mm	740mm	140kg