

## Quick Reference Guide

# Standard Shower Commode Range

Standard tilt in space mobile shower commode (KA11507)



Standard self-propel mobile shower commode (KA120SR)



Standard transit mobile shower commode (KA114S)



## Knowledge, Experience and Specialist Consultation

Years of knowledge, experience, and consultation with specialists and consumers, goes into the design and construction of KCare's premium range of shower commodes to deliver superior;

Support

Comfort

Safety

Functionality

The following quick reference guide (see over) is designed to highlight standard features, key differences and various options across 8 model ranges.

# Options and Accessories

KCare Healthcare Solutions stocks a comprehensive range of options and accessories for its shower commodes. Options can be fitted upon ordering a new commode or purchased at a later date\*. If our standard range of options do not suit your requirements, please contact a KCare Healthcare Solutions dealer for a custom-built solution.

- Contoured headrest
- Pommel splash guard
- Adjustable angle footplates
- Strapped back
- Fixed curved padded back
- Armrest pads
- Rotating safety arms
- Gutter arms
- Neoprene footplate covers
- Vinyl padded individual footplates
- Footrest extensions
- Pan
- Pan carrier
- Calf support pads
- Padded vinyl calf strap
- Stump support
- Elevating leg-rest
- Stainless steel braking castors
- Low profile swing up arms
- Tension adjustable backrest
- Wheelchair style backrest
- Tie on padded back
- Pelvic strap\*

## Quick Reference Guide

### Standard Features

- Independently tested to conform to AS/NZS 3973:2009
- Australian owned, designed and distributed
- Multiple seat widths available
- Manufactured from durable stainless steel
- High quality 125mm diameter braking castors
- Self-propel models standard with 22" wheels
- Variety of polyurethane, handmade and pressure relieving seats available (custom seats available on request)
- All models highly customisable



	Model number							
Features	KA112 Range	KA113 Range	KA114 Range	KA115 Range	KA116 Range	KA120 Range	KA123 Range	KA130 Range
<b>Total number of models in range</b>	6	2	7	6	1	10	2	3
<b>Weight Capacity (kg)</b>	200 - 300	200	120 - 300	200	200	120 - 300	200	200
<b>Paediatric Option</b>				•		•		
<b>Bariatric Option</b>	•		•	•		•		
<b>Operating Style</b>								
<b>Attendant Propelled</b>	•	•	•	•	•			
<b>Self Propelled</b>						•	•	•
<b>Tilt In Space</b>				•				
<b>Foot Rest Type</b>								
<b>Swingaway</b>			•	•	•	•		•
<b>Slide out – weight bearing</b>		•					•	
<b>Slide out – non weight bearing</b>	•							
<b>Arm Rest Type</b>								
<b>Loop Swing Up</b>	•	•	•	•	•	•	•	•
<b>Removable</b>	•		•	•		•		
<b>Back Rest Type</b>								
<b>PVC Strap</b>	•	•	•	•	•	•	•	•
<b>Padded</b>	•		•	•		•		
<b>Warranty Period</b>	12 months with 36 months structural and 5 year anti-corrosion							

\*Some options and accessories will impact other features and compatibility is dependent on the base commode model.

\*The pelvic strap is designed for positioning purposes only and is not designed to be used as a safety belt.

Please contact [customerservice@kcare.com.au](mailto:customerservice@kcare.com.au) to find out more information on specific models in each range.

[www.kcare.com.au](http://www.kcare.com.au)