



# TILT IN SPACE MOBILE SHOWER COMMUNE WITH SWINGAWAY FOOTREST - 445MM

Item No. KA11507

Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind . Extra high PVC strap backrest for improved comfort and safety . 10-35° tilt angle . Integrated push-handle for manoeuvrability . Swing-up padded arms for unrestricted side transfer . Swingaway and



## Attendant Tilt In Space Shower Commode



### DESCRIPTION

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind .
- Extra high PVC strap backrest for improved comfort and safety .
- 10-35° tilt angle .
- Integrated push-handle for manoeuvrability .
- Swing-up padded arms for unrestricted side transfer .
- Swingaway and removable height adjustable footrest with contoured padded calf supports .
- 4x 125mm total locking castors .
- Designed to be used with a K Care patented removable seats .
- Available in 400, 445mm and 500mm seat widths .
- Padded seat not included .

### SPECIFICATIONS

Overall Width	580 mm
---------------	--------

---

Overall Depth	1070 mm
---------------	---------

---

Overall Height	1453 mm
Seat Width	445 mm
Seat Depth	430 mm
Floor To Seat Height	540 mm
Armrest Type	Loop
Armrest Height	365 mm
Backrest Height	855 mm
Footrest Type	Swing Away
Product Weight	21 kg
Max User Weight	200 kg
Brand	K Care