



Everyday Living Aids Catalogue

Contents

Breezy Everyday Aqua Elite	
Shower Commode Chair, Self Propelled	5
Shower Commode Chair, Transit	5
Shower Commode Chair, Tilt N Space	6
Breezy Everyday	
Bathroom	
Bath Bench and Shower Boards	7
Bath Stools and Shower Chairs	8
Bathtub Bars	12
Transfer Benches	13
Toileting	
Raised Toilet Seats	15
Commodes	18
Urinals	20
Mobility	
Canes and Walking Sticks	21
Canes, Quad	25
Crutches, Elbow	26
Crutches, Under Arm	32
Walking Frames, Folding	36
Walking Frames, Non Folding	39
Miscellaneous	43 - 45
Spare Parts / Accessories	46 - 47

To order please contact:

Welcome to the Breezy Everyday Catalogue!

At Sunrise Medical, we pride ourselves on our commitment to supply high quality products at affordable prices.

We are pleased to be able to introduce our new Breezy Everyday Living Aids range. These practical products have been carefully designed and selected for user satisfaction.

We have a friendly and knowledgeable Customer Service Team who are there to help with any queries or problems you may have. We also have a team of Area Sales Managers who are always available for additional information, product demonstrations or quotations.

For further information, or to discuss your needs in more detail, please feel free to contact your Area Sales Manager or Customer Service Associate on (02) 9678 6600 between 9am and 5pm Monday to Friday or email enquiries@sunrisemedical.com.au.

Sunrise Medical thanks you for your continued support and looks forward to serving you in the near future.

To place your order:

E-mail: orders@sunrisemedical.com.au

To view the everyday living aids range online and to download this brochure, please visit www.sunrisemedical.com.au







Shower Commode Chair, 45cm, Hybrid

- Perfect for the everyday user or for use with an attendant
- Hybrid Self Propelled or Transit Wheels
- Aluminium powder coating
- Flip back armrest with PU cover
- Detachable swing away footrest with flip up footplates
- PU back & seat with pan and cover
- Front braked casters: 5" PU tyre
- Rear caster wheels: 6" PU tyre
- Rear wheel option: 24" x 1 3/8 PU quick release, aluminium hand rim with brake
- Attendant controlled Hand brake
- Foldable 'push to lock' wheel locks
- Intimate seat
- AU STD approved

Product Code	UOM	Specifications	Max User Weight	Warranty
BE5500 -EA	EA	Seat width: 450mm Seat depth: 420mm Width adjustable (20" extension kit available) Chair weight: 14kg (with 6" rear tyre) Stainless steel screws Fits over a toilet - std toilet height = 15.5" Commode clearance height = 17.5"	136kgs	1 year







Shower Commode Chair, Tilt N Space, Electric

- PU intimate seat with pan & cover, for better hygiene
- Flip backward armrest with PU pad
- Removable swing away hangers and flip up footplates
- Tall PU backrest
- Includes adjustable and removable headrest and laterals
- Front castors: 5" PU tyre with brake
- Rear wheels: 5" PU tyre with brake
- Stainless steel screws
- Quickie style anti-tips
- AU STD approved
- ARTG No: 224739

Product Code	UOM	Specifications	Max User Weight	Warranty
BE5502- EA	EA	Chair weight: 25.8kg (including charger 0.8kg Double electric actuators Charger: 24V / 1000mA Battery: 2 x 1.3AH x 12V Seat width: 460mm Seat depth: 425mm Tilt Range: 5° forward to 36.5° backwards Height Adjustable: 525mm to 625mm Fits over a toilet - std toilet height = 15.5" Commode clearance height = 17.5"	136kg	1 year



Shower Commode Chair, Tilt N Space, Gas

- PU intimate seat with pan & cover, for better hygiene
- Flip backward armrest with PU pad
- Removable swing away hangers and flip up footplates
- Tall PU backrest
- · Includes adjustable and removable headrest and laterals
- Front castors: 5" PU tyre with brake
- Rear wheels: 5" PU tyre with brake
- · Stainless steel screws
- Quickie style anti-tips
- Round push handle
- AU STD approved
- ARTG No: 224739

Produc Code		Specifications	Max User Weight	Warranty
BE550 EA	3- EA	2 x Stainless Steel gas springs - 300N (use for 40-80kg) 500N (use for 80-120kg) Seat width: 460mm Seat depth: 425mm Chair weight: 25kg Tilt Range: 36.5° backwards Height Adjustable: 525-625mm Fits over a toilet - std toilet height = 15.5" Commode clearance height = 17.5"	136kg	1 year






Bath Board Adjustable, PVC

- A strong shower board constructed of reinforced PVC for maximum durability
- Secured into position with a self-locating fixing system
- The slats are tapered to encourage rapid drainage
- ARTG No: 203577

Product Code	UOM	Specifications	Max User Weight	Warranty
BE10167	EA	Length: 660 mm Width: 340 mm Length Adjustment: 90 mm at each end Product Weight: 3.2kg each	190 kg	1 year



Bath Board Adjustable, Plastic

- For ease of transfer and added security, this moulded plastic bath board is supplied complete with handle
- Providing a comfortable contoured seating area for sitting whilst bathing, it has well spaced holes for improved drainage
- AU STD approved
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE68801	EA	Length: 675 mm Width: 260 mm Length Adjustment: 410 - 620 mm Product Weight: 2.42kg each	150 kg	1 year





Bariatric Bath Stool, Detach Back

- Suction-cup feet secure bench in bathtub
- Dual cross brace for extra stability
- · High-strength, aluminium frame to prevent rusting
- · Riveted for added strength
- Comes completely assembled
- Height Adjustable
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE68024HD	EA	Height Adjustment: 419 - 521 mm Seat Width: 500 mm Product Weight: 3.31kg each	225 kg	1 year



Bath Stool With Fixed Back

- High density polyethylene back and seat with 19" wide and 12" deep for maximum patient comfort and stability
- Durable aluminium frame is corrosion resistant. Angled legs distribute weight over a large area to increase patient security
- Suction-cup tips on all four legs provide extra stability and safety
- Legs quick-release design for easy transport and storage
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE68030S	EA	Seat Height: 350 - 525 mm Seat Width: 505 mm Seat Depth: 305 mm Product Weight: 3.21kg each	113 kg	1 year





Swivel Shower Chair, Powder Coated

- This shower chair has a base-frame which fits over the bath, with easy-to-clean plastic
- Sturdy
- Swivels 360°
- Secure with locking lever
- ARTG No: 202562

Prod	uct Code	UOM	Specifications	Max User Weight	Warranty
BE	68039	EA	Backrest Height: 356 mm Seat Width: 415 mm Seat Depth: 410 mm Distance between armrests: 490 mm Product Weight: 6kg each	113 kg	1 year



Shower Stool Without Back

Open front plastic moulded seat

ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE68054	EA	Seat Width: 410 mm Seat Depth: 400 mm Height Adjustment: 432 - 508 mm Arms Width: 470 mm Arm Height: 220 mm	113 kg	1 year





Shower Chair with Adjustable Height and Armrests

- The plastic seat is curved to give comfort and has drainage holes for excess water
- AU STD approved
- ARTG No: 202562

	Product Code	UOM	Specifications	Max User Weight	Warranty
BE	E68116H	EA	Seat Width: 420 mm Seat Depth: 400 mm Width between arms: 540 mm Height Adjustment: 457 - 635 mm	150 kg	1 year





Fixed Height Shower Chair with Wheels, Arms & Hangers

- The plastic seat is curved to give comfort and has drainage holes for excess water
- Detachable legs
- 2 Braked Casters
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3656	EA	Seat Width: 420 mm Seat Depth: 400 mm Seat Height: 490 mm Arm Width: 540mm	113 kg	1 year





Bathtub Bar, White

- White enamel-coated 1" steel bar provides firm handle for reliable patient assistance and increased safety during transfers
- Adjustable clamp attaches to tub walls up to 7" wide
- Unit is 6" wide and extends 14" above the tub wall
- Equipped with rubber pads to protect walls
- Rubber pads for security and protection
- ARTG No: 160747

Product Code	UOM	Specifications	Max User Weight	Warranty
BE67050	EA	Width: 6" Height: 14"	190 kg	1 year





Unpadded Bariatric Transfer Bench, Assembled

- Textured, moulded backrest and seating surface offer a non-slip surface with holes for water drainage
- Aluminium frame is height adjustable in 1/2" increments
- Backrest adjusts without tools for right or left facing tubs
- Side arm and seat handle provide extra stability and leverage when transferring
- Safely provides assistance to those who have difficulty stepping over bathtub walls
- Backrest and armrest installed with detent buttons for tool free assembly
- Has stainless steel spring pins
- ARTG No: 203577

Product Code	UOM	Specifications	Max User Weight	Warranty
BE98310B	EA	Seat Width: 670 mm Height Adjustment: 390 - 520 mm	227 kg	1 year



Unpadded Transfer Bench, Assembled

- Suction cup foot pieces secure the bench to tub floor
- Lightweight, anodised aluminium frame is height adjustable in 1/2" increments to accommodate most tubs
- · Backrest adjusts without tools for right or left facing tubs
- Built-in dish and hand-held shower slot for user convenience
- Has stainless steel spring pins
- Extension legs (80 mm) available (see BE3659, BE3660, BE3661)
- ARTG No: 203577

Product Code	UOM	Specifications	Max User Weight	Warranty
BE98309	EA	Seat Width: 800 mm Seat Depth: 370 mm Height Adjustment: 430 - 550 mm Product Weight: 2kg each	136 kg	1 year





Transfer Bench, Deluxe Padded, Assembled

- Extra strong, watertight heat sealed vinyl padding on seat and back is comfortable, easy to clean and durable
- · Super hard anodised aluminium frame tubing
- Wide base for extra stability
- Suction cup tips for two legs of the bath bench provide extra stability and safety
- · Has stainless steel spring pins
- · Detachable backrest and armrest
- ARTG No: 160747

Product Code	UOM	Specifications	Max User Weight	Warranty
BE68011P	EA	Seat Width: 610 mm Seat Depth: 410 mm Height Adjustment: 445 - 546 mm Product Weight: 5.2kg each	113 kg	1 year



Transfer Bench, Padded Seat & Backrest, Assembled

- Vinyl covered foam backrest & seating surface offer a highly resilient cushioned surface for comfort.
- Lightweight, anodised aluminium frame is height adjustable in 1/2" increments.
- Suction cup foot pieces secure the bench to the tub floor.
- · Backrest adjusts for right or left facing tubs
- Side arm & seat handle provide extra stability & leverage when transferring
- Has stainless steel spring pins
- ARTG No: 203577

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3645	EA	Seat Width: 685 mm Height Adjustment: 420 - 552 mm Product Weight: 4.32kg each	136 kg	1 year





Raised Toilet Seat Without Lid 2"

- · Lightweight, all plastic construction, for strength and durability
- The seat has a gentle contoured surface for extra comfort
- Resistant to stains and odours, it is fully sealed and easy to clean
- The large front and rear cut-outs allow for easier personal cleansing
- It has no front bracket or front turning knob, making fitting easier and reducing the risk of snagging skin or clothes
- AU STD approved
- ARTG No: 202562

Produc	ct Code	UOM	Specifications	Max User Weight	Warranty
BE6	7212	EA	Length: 407 mm Width: 369 mm Height: 51 mm	190 kg	1 year



Raised Toilet Seat Without Lid 4"

- · Lightweight, all plastic construction, for strength and durability
- The seat has a gentle contoured surface for extra comfort
- Resistant to stains and odours, it is fully sealed and easy to clean
- The large front and rear cut-outs allow for easier personal cleansing
- It has no front bracket or front turning knob, making fitting easier and reducing the risk of snagging skin or clothes
- AU STD approved
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE67213	EA	Length: 407 mm Width: 369 mm Height: 102 mm	190 kg	1 year





Raised Toilet Seat Without Lid 6"

- Lightweight, all plastic construction, for strength and durability
- The seat has a gentle contoured surface for extra comfort
- · Resistant to stains and odours, it is fully sealed and easy to clean
- The large front and rear cut-outs allow for easier personal cleansing
- It has no front bracket or front turning knob, making fitting easier and reducing the risk of snagging skin or clothes
- AU STD approved
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE67214	EA	Length: 407 mm Width: 369 mm Height: 152 mm	190 kg	1 year



Raised Toilet Seat With Lid 2"

- Lightweight with lid, all plastic construction, for strength and durability
- The seat has a gentle contoured surface for extra comfort
- Resistant to stains and odours, it is fully sealed and easy to clean
- The large front and rear cut-outs allow for easier personal cleansing
- It has no front bracket or front turning knob, making fitting easier and reducing the risk of snagging skin or clothes
- AU STD approved
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE67212L	EA	Length: 407 mm Width: 369 mm Height: 51 mm	190 kg	1 year





Raised Toilet Seat With Lid 4"

- Lightweight with lid, all plastic construction, for strength and durability
- · The seat has a gentle contoured surface for extra comfort
- Resistant to stains and odours, it is fully sealed and easy to clean
- The large front and rear cut-outs allow for easier personal cleansing
- It has no front bracket or front turning knob, making fitting easier and reducing the risk of snagging skin or clothes
- AU STD approved
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE67213L	EA	Length: 407 mm Width: 369 mm Height: 102 mm	190 kg	1 year



Raised Toilet Seat With Lid 6"

- · Lightweight with lid, all plastic construction, for strength and durability
- The seat has a gentle contoured surface for extra comfort
- Resistant to stains and odours, it is fully sealed and easy to clean
- The large front and rear cut-outs allow for easier personal cleansing
- It has no front bracket or front turning knob, making fitting easier and reducing the risk of snagging skin or clothes
- AU STD approved
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE67214L	EA	Length: 407 mm Width: 369 mm Height: 152 mm	190 kg	1 year





3 in 1 Adjustable Bariatric Commode

- · Back rail removes quickly and easily
- Features a durable, plastic seat with lid
- Plastic armrests for additional patient comfort
- Includes 10-quart capacity, commode bucket with cover and splash shield
- Non-skid, non-marring, rubber tips provide extra safety
- Adjustable legs allow seat to be raised or lowered in 1" increments to sit from 496 mm - 597 mm (19.5" - 23.5") bariatric
- Fits over a toilet
- ARTG No: 202562

Product Code	e UOM	Specifications	Max User Weight	Warranty
BE68111	EA	Bariatric: 496 - 597 mm Width between arms: 590mm Height Adjustment: 496 - 597 mm, in 1" increments	225 kg	1 year



Extra Wide Bariatric Commode, With Adjustable Padded Back

- All steel frame with Versa Guard coating is durable and warm to the touch
- Extra-side platform seat maximises stability and is comfortable for long periods of sitting
- · Comes complete with pail and lid.
- Fits over a toilet
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE68115XW	EA	Height Adjustment: 496 - 597 mm Width between arms: 580 mm Minimum Height: 430 mm - 630 mm	225 kg	1 year





Deluxe 3 in 1, Heavy Duty Adjustable Commode

- Height Adjustable
- Removable solid plastic toilet seat and cover
- Full length swing- away rubber armrests provide comfort
- Powder coated steel
- Fits over a toilet
- Swing-away arms
- Material: Steel
- ARTG No: 202562

Product Code	UOM	Specifications	Max User Weight	Warranty
BE68680DA	EA	Armrest Height: 6.5" (above seat) Width: 550 mm Width between arms: 460 mm Height from arms to floor: 673 mm - 775 mm Seat Height Adjustment: 485 - 590 mm	225 kg	1 year





Male Urinal With Cover

• Easy grip open handle

Hangs on bed rail

• ARTG No: 202562

Product Code	UOM	Capacity	Specifications	Warranty
BE28101	EA	1000 cc	Material:Translucent HDPE body with gold cap Product Weight: 135g each	1 year



Urinal, Female

Translucent HDPE

• ARTG No: 202562

Product Code	UOM	Capacity	Specifications	Warranty
BE28102	EA	1000 cc	Material: Translucent HDPE Product Weight: 140g each	1 year



Urinal, Female, Fractured Bed Pan

• Polypropylene

Autoclavable

• ARTG No: 202562

Product Code	UOM	Dimensions	Specifications	Warranty
BE28202	EA	330 x 230 x 80mm	Material: Polypropylene Product Weight: 245g each	1 year



Urinal, Female, Fractured Bed Pan, Large

Polypropylene

Autoclavable

• ARTG No: 202562

Product Code	UOM	Dimensions	Specifications	Warranty
BE3657	EA	437 x 305 x 80mm	Material: Polypropylene Product Weight: 502g each	1 year





Adjustable Fischer Stick, Black, Left Hand, Medium

- A unique, biometrically shaped handle designed to evenly space pressure across the palm of the hand
- Useful for arthritics or those with poor grip
- AU STD approved
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1655	EA	Left hand: 670 - 925 mm Ferrule: 3/4", Black (BE3653) Product Weight: 400g each	127 kg	1 year



Adjustable Fischer Stick, Black, Right Hand, Medium

- A unique, biometrically shaped handle designed to evenly space pressure across the palm of the hand
- Useful for arthritics or those with poor grip
- AU STD approved
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1657	EA	Right hand: 670 - 925 mm Ferrule: 3/4", Black (BE3653)Product Weight: 400g each	127 kg	1 year



Cane With Plastic Handle

- Height adjustable support stick fitted with ultra-safe horseshoe clips for quick and easy adjustability
- AU STD approved
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1665P	EA	Height Adjustment: 715 - 970 mm Ferrule: 3/4", Black (BE3653 Product Weight: 365g each	100 kg	1 year





Cane Foam Swan Neck With Strap

- Height adjustable support stick fitted with ultra-safe horseshoe clips for quick and easy adjustability
- AU STD approved
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1674	EA	Height Adjustment: 730 - 985 mm Ferrule: 3/4", Black (BE3653) Product Weight: 365g each	100 kg	1 year



Aluminium Cane, Swan Neck With PVC Handle

- A strong, lightweight support stick which is height adjustable
- Swan neck which helps to transmit the weight of the user vertically through the body of the stick
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1711	EA	Height Adjustment: 710 - 965 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 550g each	127 kg	1 year



Aluminium Cane, Straight Neck With PVC Handle, Medium

- Straight neck PVC handle
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3001	EA	Height Adjustment: 710 - 965 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 550g each	127 kg	1 year





Aluminium Cane Straight Neck, Medium

Straight neck, plastic handle ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3004	EA	Height Adjustment: 700 - 960 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 450g each	127 kg	1 year



Aluminium Cane, Straight Neck With Plastic Handle, Large

Straight neck, plastic handle ARTG No: 160745

	oduct Code	UOM	Specifications	Max User Weight	Warranty
BE	3006	EA	Height Adjustment: 845 - 1100 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 450g each	127 kg	1 year






Bariatric Cane, Offset Handle, Black

- Solid brass, push button adjustment for extra stability
- Handle height adjusts in 1" increments from 737 966mm (29" 38")
- Secured with a locking ring for added safety

Ferrule: 22 mm ARTG No: 203275

Produc Code		Specifications	Colour	Max User Weight	Warranty
BE710	66 EA	Height Adjustment: 737 - 966 mm Ferrule: 22.1 mm	Black	225 kg	1 year



Aluminium Cane, Standard Handle, Adjustable, Chrome

- Contoured handles with foam hand grips provide a comfortable, secure grip
- Attractive finishes resist marring and scratching, while providing a non-institutional look
- Handle height adjusts in 1" increments
- Secured with a locking ring for added safety
- Vinyl handgrip Ferrule: 19 mm
- ARTG No: 203275

Product Code	UOM	Specifications	Colour	Max User Weight	Warranty
BE72092	EA	Height Adjustment: 762 - 991 mm Ferrule: 19 mm	Chrome	113 kg	1 year



Folding Cane, Contour Handle, Strap, Black

- Quick and simple to use, this walking stick opens in seconds
- Lightweight yet strong and stable, these high quality folding walking sticks conveniently fold away into a plastic wallet to slip into a bag, briefcase etc
- Height adjustable Ferrule: 19 mm ARTG No: 160745

Product Code	UOM	Specifications	Colour	Max User Weight	Warranty
BE9889	EA	Height Adjustment: 787 - 900 mm Ferrule: 3/4", Black (BE3653)	Black	100 kg	1 year





Aluminium Quad Cane, Black, Large Base, Adult

- Available in black
- Engineered and tested to ensure performance and quality
- Soft foam handgrip enhances comfort
- Locking silencer provides rattle free use
- · Solid brass detent button resists wear and breakage
- Height adjusts easily
- ARTG No: 203275

Product Code	UOM	Specifications	Colour	Max User Weight	Warranty
BE30341	EA	User Height: 1470 - 1930 mm Base Size: 200 x 300mm Ferrule: 16 mm (BE3670)	Black	113 kg	1 year



Aluminium Quad Cane, Black, Small Base, Adult

- Available in black
- Engineered and tested to ensure performance and quality
- Soft foam handgrip enhances comfort
- Locking silencer provides rattle free use
- Solid brass detent button resists wear and breakage
- Height adjusts easily
- ARTG No: 203275

Product Code	UOM	Specifications	Colour	Max User Weight	Warranty
BE30346	6 EA	User Height: 1170 - 1630 mm Base Size: 150 x 200mm Ferrule: 13 mm (BE3650)	Black	113 kg	1 year



Aluminium Quad Cane, Deluxe, Large Base, Adult

- Available in black
- Engineered and tested to ensure performance and quality
- Soft foam handgrip enhances comfort
- · Locking silencer provides rattle free use
- Solid brass detent button resists wear and breakage
- · Height adjusts easily
- Ferrule: 16 mm
- ARTG No:160745

Product Code	UOM	Specifications	Colour	Max User Weight	Warranty
BE2620BL	EA	Height Adjustment: 775 - 1025 mm User Height: 1550 - 2050mm Base size: 200 x 300 mm Ferrule: 5/8", Grey (BE3651) Product Weight: 1.5kg each	Black	113 kg	1 year





Children's Elbow Crutch, 4 - 7 Years

A cheerful buttercup yellow elbow crutch designed for children.

- Made from strong yet lightweight aluminium
- · Height adjustable
- Latex free (including ferrule)
- ARTG No:160745

Product Code	UOM	Specifications	Colour	Max User Weight	Warranty
BE3640	Pair	Handgrip Height: 406 - 559 mm Ferrule: 3/4", Grey (BE3652)	Yellow	100 kg	1 year



Children's Elbow Crutch, 6 - 10 Years

A cheerful yellow elbow crutch designed for children.

- Made from strong yet lightweight aluminium
- · Height adjustable
- Latex free (including ferrule)
- ARTG No:160745

Product Code	UOM	Specifications	Colour	Max User Weight	Warranty
BE1070	Pair	Handgrip Height: 457 - 656 mm Ferrule: 3/4", Grey (BE3652)	Yellow	100 kg	1 year





Arthritic Elbow Crutch

For arthritics or those with impaired grip.

- Adjustable height
- A padded forearm trough distributes the user's weight over a large support area
- The adjustable telescopic handle rotates through 360° and locks to suit individual forearm length
- Latex free (including ferrules)
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1300	Pair	Trough Height: 980 - 1235 mm Ferrule: 1", Grey (BE3067) Product Weight: 950g each	127 kg	1 year



Plastic Handle Elbow Crutch Closed Cuff, Double Adjustable, Medium

A strong lightweight double adjustable crutch with a shaped hand grip and arm support.

- Height adjustable aluminium crutches
- The shape of the handle on this elbow crutch has been designed to spread the pressure more evenly over the hand. The triangular shape and angle of the handle reflect the natural shape of the hand
- Double adjustable crutches are adjustable between cuff and handle and handle and ground
- Latex free (including ferrule)
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1263	Pair	Handgrip Height: 530 - 785 mm Cuff to Handgrip: 200 - 275 mm Cuff Diameter: 95 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 600g each	180 kg	1 year





Plastic Handle Elbow Crutch Closed Cuff, Double Adjustable, Adult

A strong lightweight double adjustable crutch with a shaped hand grip and arm support.

- · Height adjustable aluminium crutches
- The shape of the handle on this elbow crutch has been designed to spread the pressure more evenly over the hand. The triangular shape and angle of the handle reflect the natural shape of the hand
- Double adjustable crutches are adjustable between cuff and handle and handle and ground
- Latex free (including ferrule)
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1265	Pair	Handgrip Height: 635 - 910 mm Cuff to Handgrip: 200 - 295 mm Cuff Diameter: 95 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 600g each	180 kg	1 year



Plastic Handle Elbow Crutch Closed Cuff, Double Adjustable, Extra Long

A strong lightweight double adjustable crutch with a shaped hand grip and arm support.

- · Height adjustable aluminium crutches
- The shape of the handle on this elbow crutch has been designed to spread the pressure more evenly over the hand. The triangular shape and angle of the handle reflect the natural shape of the hand
- Double adjustable crutches are adjustable between cuff and handle and handle and ground
- Latex free (including ferrule)
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1267	Pair	Handgrip Height: 820 - 1075 mm Cuff to Handgrip: 200 - 295 mm Cuff Diameter: 95 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 600g each	180 kg	1 year





Plastic Handle Elbow Crutch Closed Cuff, Single Adjustable, Adult

A strong lightweight single adjustable crutch with a shaped hand grip and arm support.

- Height adjustable aluminium crutches
- The shape of the handle on this elbow crutch has been designed to spread the pressure more evenly over the hand. The triangular shape and angle of the handle reflect the natural shape of the hand
- Single adjustable crutches are adjustable between handle and ground
- Latex free (including ferrule)
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3098	Pair	Height Adjustment: 27" - 37" Handgrip Height: 660 - 915 mm Cuff to Handgrip: 240 mm Cuff Diameter: 95 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 600g each	180 kg	1 year



Plastic Handle Elbow Crutch Closed Cuff, Single Adjustable, Adult Tall

A strong lightweight double adjustable crutch with a shaped hand grip and arm support.

- Height adjustable aluminium crutches
- The shape of the handle on this elbow crutch has been designed to spread the pressure more evenly over the hand, and the triangular shape and angle of the handle reflect the natural shape of the hand
- Single adjustable crutches are adjustable between handle and ground
- Latex free (including ferrule)
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3099	Pair	Height Adjustment: 33" - 43" Handgrip Height: 838 - 1092 mm Cuff to Handgrip: 240 mm Cuff Diameter: 95 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 600g	180 kg	1 year





Comfy Handle Elbow Crutch Closed Cuff, Double Adjustable, Medium

- This crutch offers a unique patented handle which offers pressure distribution characteristics greatly superior to any traditional crutch
- Latex free (including ferrule)
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1351	Pair	Handgrip Height: 580 - 835 mm Cuff to Handgrip Height: 220 - 275 mm Cuff Diameter: 95 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 715g each	127 kg	1 year



Comfy Handle Elbow Crutch Closed Cuff, Double Adjustable, Adult

- This crutch offers a unique patented handle which offers pressure distribution characteristics greatly superior to any traditional crutch
- Latex free (including ferrule)
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1353	Pair	Handgrip Height: 680 - 935 mm Cuff to Handgrip Height: 220 - 275 mm Cuff Diameter: 95 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 715g each	127 kg	1 year





Permanent User Elbow Crutch, PVC Handle

- This extra strong crutch made from reinforced aluminium is designed for those users with longer term needs. The absence of adjustment holes makes the crutch "silent" in use. The crutch is supplied in one standard length, and may be cut to the required length. It comes with no ferrule- order separately (BE3067).
- Latex free (including ferrules)
- Cut to length
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE1043	Pair	Handle to Floor Height: 960 mm Cuff to Handle: 229 mm Cuff Diameter: 100 mm Ferrule: 1", Grey (BE3067) Product Weight: 800g each	190 kg	1 year



Permanent User Elbow Crutch, Comfy Handle

- This extra strong crutch made from reinforced aluminium is designed for those users with longer term needs. The absence of adjustment holes makes the crutch "silent" in use. The crutch is supplied in one standard length, and may be cut to the required length. It comes with no ferrule- order separately (BE3067).
- Latex free (including ferrules)
- Cut to length
- ARTG No:160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3042	Pair	Cuff Diameter: 100 mm Handle to Floor Height: 960 mm Cuff to Handle: 229 mm Ferrule: 1", Grey (BE3067) Product Weight: 800g each	190 kg	1 year





Children's Underarm Crutch, 4 - 7 Years


- · Easy to use
- Latex free (including ferrules)
- ARTG No:160745



Product Code	UOM	Specifications	Max User Weight	Warranty
BE3641	Pair	Handgrip to Armrest: 178 - 292 mm Floor to Armrest: 787 - 940 mm Ferrule: 7/8", Grey (BE3066) Product Weight: 1.42kg / pair	136 kg	1 year

3 in 1 crutch

- · Anodised aluminium silver tube (Nature Silver)
- One size fits 4'7" 6'7" individuals
- · Foam arm pad and handgrip enhances comfort
- Non-skid rubber tip
- Durable and lightweight
- Easy to use
- · Spare parts include cushions, hand grips and tips
- Latex free
- ARTG No:160745







Red Dot Crutch, Aluminium, Tall

- Clip-lock adjusted handle and push-button adjusted height allows for 10 second adjustments
- 1" adjustment increments on both handle and height, provide a precise fit which reduces muscle fatigue
- High density latex free rubber cushions and hand grips provide long wear and greater comfort
- Latex free
- AU STD approved
- ARTG No: 203275

Product Code	UOM	Specifications	Max User Weight	Warranty
BE90-214-8	Pair	Height Adjustment: 1320 - 1520 mm User Height: 1780 - 1980 mm Ferrule: 22.6 mm Crutch Cushion: BE3646 Crutch Handgrip: BE3654	136 kg	1 year



Red Dot Crutch, Aluminium, Adult

- Clip-lock adjusted handle and push-button adjusted height allows for 10 second adjustments
- 1" adjustment increments on both handle and height, provide a precise fit which reduces muscle fatigue
- High density latex free rubber cushions and hand grips provide long wear and greater comfort
- Latex free
- · AU STD approved
- ARTG No: 203275

Product Code	UOM	Specifications	Max User Weight	Warranty
BE91-214-8	Pair	Height Adjustment: 1120 - 1320 mm User Height: 1550 - 1750 mm Ferrule: 22.6 mm Crutch Cushion: BE3646 Crutch Handgrip: BE3654	136 kg	1 year





Red Dot Crutch, Aluminium, Youth

- Clip-lock adjusted handle and push-button adjusted height allows for 10 second adjustments
- 1" adjustment increments on both handle and height, provide a precise fit which reduces muscle fatigue
- High density latex free rubber cushions and hand grips provide long wear and greater comfort
- Latex free
- AU STDS approved
- ARTG No: 203275

Product Code	UOM	Specifications	Max User Weight	Warranty
BE92-214-8	Pair	Height Adjustment: 940 - 1140 mm User Height: 1320 - 1570 mm Ferrule: 22.6 mm Crutch Cushion: BE3646 Crutch Handgrip: BE3654	136 kg	1 year



Bariatric Crutch, Steel, Adult

- Heavy-duty steel construction provides durability and supports up to 249kg
- Push-button adjusted foot pieces with fitting scale allow for accurate, fast fitting
- 1" adjustment increments on both handle and height, provide a precise fit which reduces muscle fatigue
- Latex free
- ARTG No: 203275

Product Code	UOM	Specifications	Max User Weight	Warranty
BE61314B	Pair	Height Adjustment: 1575 - 1778 mm Ferrule: 22.2 mm Crutch Cushion: BE3647 Crutch Handgrip: BE3649	249 kg	1 year




Crutches, Under Arm (Axilla)



Bariatric Crutch, Steel, Tall

- Heavy-duty steel construction provides durability and supports up to 249kg
- Push-button adjusted foot pieces with fitting scale allow for accurate, fast fitting
- 1" adjustment increments on both handle and height, provide a precise fit which reduces muscle fatigue
- Latex free
- ARTG No: 203275

Product Code	UOM	Specifications	Max User Weight	Warranty
BE60314B	Pair	Height Adjustment: 1778 - 1981 mm Ferrule: 22.2 mm Crutch Cushion: BE3647 Crutch Handgrip: BE3649	249 kg	1 year





Folding, Small, W/O Wheels

- A simple folding frame which does not use buttons, catches or other tricky mechanisms. To fold simply pull the ball towards you
- When folded, the frame fits easily into small storage spaces
- Latex-free
- · Lightweight aluminium
- No wheels
- Height Adjustable
- · Frame comes already assembled
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3621A	EA	Maximum Width: 630 mm Maximum Depth: 650 mm Handgrip Height: 711 - 787mm Ferrule: 1" Grey (BE3067) Product Weight: 1.8kg each	127 kg	1 year



Folding, Medium, W/O Wheels

- A simple folding frame which does not use buttons, catches or other tricky mechanisms. To fold simply pull the ball towards you
- When folded, the frame fits easily into small storage spaces
- Latex-free
- Lightweight aluminium
- No wheels
- Height Adjustable
- · Frame comes already assembled
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3622A	EA	Maximum Width: 635 mm Maximum Depth: 555 mm Handgrip Height: 787 - 863mm Ferrule: 1" Grey (BE3067) Product Weight: 2kg each	127 kg	1 year



Folding, Large, W/O Wheels

- A simple folding frame which does not use buttons, catches or other tricky mechanisms. To fold simply pull the ball towards you
- When folded ,the frame fits easily into small storage spaces
- Latex-free
- · Lightweight aluminium
- No wheels
- Height Adjustable
- · Frame comes already assembled
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3623A	EA	Maximum Width: 635 mm Maximum Depth: 654 mm Handgrip Height: 838 - 914 mm Ferrule: 1" Grey (BE3067) Product Weight: 2.2kg each	127 kg	1 year



Folding, Medium, With Wheels

- A simple folding frame which does not use buttons, catches or other tricky mechanisms - to fold simply pull the ball towards you
- When folded the frame fits easily into small storage space
- · Latex-free
- · Lightweight aluminium
- Wheels
- Height Adjustable
- · Frame comes already assembled
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3625A	EA	Maximum Width: 640 mm Maximum Depth: 570 mm Height Adjustable: 787 - 863 mm Ferrule: 1", Grey (BE3067) Product Weight: 2.1kg each	127 kg	1 year



Folding, Large, With Wheels

- A simple folding frame which does not use buttons, catches or other tricky mechanisms - to fold simply pull the ball towards you
- When folded the frame fits easily into small storage space
- Latex-free
- · Lightweight aluminium
- Wheels
- Height Adjustable
- · Frame comes already assembled
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3626A	EA	Maximum Width: 640 mm Maximum Depth: 570 mm Height Adjustable: 838 - 914 mm Ferrule: 1", Grey (BE3067) Product Weight: 2.3kg each	127 kg	1 year





Folding Walking Frame, Heavy Duty, W/O Wheels

- · Heavy duty walker designed with durability and style
- Extra wide frame, higher weight capacity and more durable components for increased comfort and stability
- Height Adjustable
- No Wheels
- · Frame comes already assembled
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3637	EA	Maximum Width: 740 mm Depth Open: 600 mm Depth Folded: 100 mm Height Adjustable: 810 - 990 mm Product Weight: 8.3kg each	227 kg	1 year



Folding Walking Frame, Heavy Duty, Wheels

- · Heavy duty walker designed with durability and style
- Extra wide frame, higher weight capacity and more durable components for increased comfort and stability
- Wheels
- Height Adjustable
- · Frame comes assembled
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3638	EA	Maximum Width: 740 mm Depth Open: 600 mm Depth Folded: 100 mm Height Adjustable: 810 - 990 mm Ferrule: 1-1/8" (BE3658) Product Weight: 8.3kg each	227 kg	1 year



Non Folding, Domestic Frame, Small, W/O Wheels

- High tensile extruded aluminium with profiled handgrip for user safety and comfort
- · Designed for use indoors
- No Wheels
- Height Adjustable
- · Assembly required
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3627	EA	Maximum Width: 560 mm Maximum Depth: 520 mm Height Adjustable: 665 - 740 mm User Height: 1330 - 1480 mm Ferrule: 1", Grey (BE3067) Product Weight: 1.9kg each	160 kg	1 year



Non Folding, Domestic Frame, Medium, W/O Wheels

- High tensile extruded aluminium with profiled handgrip for user safety and comfort
- · Designed for use indoors
- No Wheels
- · Height Adjustable
- · Assembly required
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3628	EA	Maximum Width: 560 mm Maximum Depth: 520 mm Height Adjustable: 765 - 840 mm User Height: 1530 - 1680 mm Ferrule: 1", Grey (BE3067) Product Weight: 1.9kg each	160 kg	1 year



Non Folding, Domestic Frame, Large, W/O Wheels

- High tensile extruded aluminium with profiled handgrip for user safety and comfort
- Designed for use indoors
- No Wheels
- Height Adjustable
- Assembly required
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3629	EA	Maximum Width: 560 mm Maximum Depth: 520 mm Height Adjustable: 865 - 940 mm User Height: 1730 - 1880 mm Ferrule: 1", Grey (BE3067) Product Weight: 1.9kg each	160 kg	1 year





Non Folding Hospital Frame, Medium, W/O Wheels

- Designed for use in hospitals and nursing homes featuring a wide and extra stable base
- No Wheels
- Height Adjustable
- · Assembly required
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3633	EA	Maximum Width: 660 mm Maximum Depth: 585 mm Height Adjustable: 765 - 840 mm User Height: 1530 - 1680 mm Ferrule: 1", Grey (BE3067) Product Weight: 1.9kg each	160 kg	1 year



Non Folding Hospital Frame, Large, W/O Wheels

- Designed for use in hospitals and nursing homes featuring a wide and extra stable base
- No Wheels
- Height Adjustable
- Assembly required
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3634	EA	Maximum Width: 660 mm Maximum Depth: 585 mm Height Adjustable: 865 - 940 mm User Height: 1730 - 1880 mm Ferrule: 1", Grey (BE3067) Product Weight: 1.9kg each	160 kg	1 year



Non Folding Hospital Frame, Medium, Wheels

- Designed for use in hospitals and nursing homes featuring a wide and extra stable base
- Wheels
- Height Adjustable
- · Assembly required
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3635	EA	Maximum Width: 660 mm Maximum Depth: 585 mm Height Adjustable: 765 - 840 mm User Height: 1530 - 1680 mm Ferrule: 1", Grey (BE3067) Product Weight: 1.9kg each	160 kg	1 year





Non Folding Hospital Frame, Large, Wheels

- Designed for use in hospitals and nursing homes featuring a wide and extra stable base
- Wheels
- Height Adjustable
- · Assembly required
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3636	EA	Maximum Width: 660 mm Maximum Depth: 585 mm Height Adjustable: 865 - 940 mm User Height: 1730 - 1880 mm Ferrule: 1", Grey (BE3067) Product Weight: 2kg each	160 kg	1 year



Non Folding, Domestic Frame, Small, Wheels

- High tensile extruded aluminium with profiled handgrip for user safety and comfort
- Designed for use indoors
- Wheels
- Height Adjustable
- · Assembly required
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3630	EA	Maximum Width: 560 mm Maximum Depth: 520 mm Height Adjustable: 665 - 740 mm User Height: 1330 - 1480 mm Ferrule: 1", Grey (BE3067) Product Weight: 2kg each	160 kg	1 year



Non Folding, Domestic Frame, Medium, Wheels

- High tensile extruded aluminium with profiled handgrip for user safety and comfort
- Designed for use indoors
- Wheels
- Height Adjustable
- · Assembly required
- ARTG No: 160745

Product Code	UOM	Specifications	Max User Weight	Warranty
BE3631	EA	Maximum Width: 560 mm Maximum Depth: 520 mm Height Adjustable: 765 - 840 mm User Height: 1580 - 1680 mm Ferrule: 1", Grey (BE3067) Product Weight: 2kg each	160 kg	1 year





Non Folding, Domestic Frame, Large, Wheels

- High tensile extruded aluminium with profiled handgrip for user safety and comfort
- Designed for use indoors
- Wheels
- Height Adjustable
- Assembly required
- ARTG No: 160745

	roduct Code	UOM	Specifications	Max User Weight	Warranty
BI	E3632	EA	Maximum Width: 560 mm Maximum Depth: 520 mm Height Adjustable: 865 - 940 mm User Height: 1730 - 1880 mm Ferrule: 1", Grey (BE3067) Product Weight: 2kg each	160 kg	1 year





Pedal Exerciser, Fixed

- Heavy duty steel frame
- Fixed
- Smooth pedal movement
- ARTG No: 219370

Product Co	ode	UOM	Specifications	Warranty
BE26001	F	EA	Colour: Blue Product Weight: 0.8kg each	1 year



Reacher, Economy

- An easy to use aid which helps the user to pick up items without bending or crouching
- ARTG No: 219370

Product Code	UOM	Specifications	Warranty
BE3642	EA	Maximum Length: 660 mm Product Weight: 0.72kg each	1 year





Long Handle Sponge

- · The long handle sponge allows for safer bathing
- ARTG No: 202562


Prod	uct Code	UOM	Specifications	Warranty
В	E3643	EA	Maximum Length: 560 mm Maximum Width: 120 mm Product Weight: 0.6kg each	1 year



Cord Style Stocking and Sock Aid

- Aids the user to put on socks, stockings or tights easily without the need to bend down.
- Simply roll the sock, stockings or tights onto the aid, slip in your foot, then use the straps to pull the clothing item up your leg
- ARTG No: 219370

Product Code	UOM	Specifications	Warranty
BE3644	EA	Maximum Length: 230 mm Maximum Width: 135 mm Product Weight: 0.8kg each	1 year



Oxygen Tank Holder

- Universal Oxygen Tank Holder
- Suited to 207, 406 and 682 size cylinders
- Will fit: Basix 2, Rubix, Breezy EC, X-series power chairs, Q2 and Q2 Lite
- ARTG No: 199075

Product Code	UOM	Specifications	Colour	Warranty
BEOXY	EA	Minimum Length: 410 mm Maximum Length: 710 mm Maximum Width: 160 mm Height Adjustable Range: 410 mm - 1080 mm Product Weight: 1.4kg each	Powder coated black	1 year



Spare Parts / Accessories



Crutch Cushion, Grey, for 3 in 1 Crutches

- Designed with thick cushioning for extra comfort
- Durable latex free rubber
- Fits BE2005
- Product code: BE3639
- UOM: Pair

Crutch Cushion, Grey, Axilla

- · Designed with thick cushioning for extra comfort
- Durable latex free rubber
- Fits BE90-214-8, BE91-214-8, BE92-214-8
- Product code: BE3646
- UOM: Pair

Crutch Cushion, Grey, Axilla, Bariatric

- Designed with thick cushioning for extra comfort
- Durable latex free rubber
- Fits BE61314B & BE60314B
- Product code: BE3647
- UOM: Pair

Crutch Hand Grip, Grey, 3 in 1

- Durable latex-free
- · Rubber provides a comfortable and sturdy grip
- Fits BE2005
- Product code: BE3648
- UOM: Pair

Crutch Hand Grip, Grey, Axilla, Bariatric

- Durable latex-free
- Rubber provides a comfortable and sturdy grip
- Fits BE61314B & BE60314B
- Product code: BE3649
- UOM: EA

Crutch Hand Grip, Red Dot

- Durable latex-free
- Rubber provides a comfortable and sturdy grip
- Fits BE90-214-8, BE91-214-8 & BE92-214-8
- Product code: BE3654
- UOM: EA

Ferrule, Rubber, Grey **Type 7/8"**

- Latex free
- High quality, moulded rubber ferrules
- Product code: BE3066
- UOM: EA

Spare Parts / Accessories



Ferrule, Rubber, Grey Type 1" (25mm)

- Latex free
- · High quality, moulded rubber ferrules
- Product code: BE3067
- UOM: EA

Ferrule, Quad Cane (Small Base), Black

Type: 13mm

- Latex-free
- Rubber tip flexes for cushioning and shock absorption
- Large surface are provides superior tracking
- High quality, moulded rubber ferrules
- Fits BE30346
- Product code: BE3650

Ferrule, Quad Cane (Large Base), Black

Type: 16mm

- Latex-free
- Rubber tip flexes for cushioning and shock absorption
- Large surface are provides superior tracking
- · High quality, moulded rubber ferrules
- Fits BE30341
- Product code: BE3670

Ferrule, Rubber, Grey Type 5/8" (16mm)

Latex-free

- Rubber tip flexes for cushioning and shock absorption
- Large surface are provides superior tracking
- · High quality, moulded rubber ferrules
- · Reinforced with metal insert for extended wear
- Product code: BE3651
- UOM: EA

Ferrule, Rubber, Grey Type 3/4" (18mm)

- Latex-free
- Rubber tip flexes for cushioning and shock absorption
- Large surface are provides superior tracking
- · High quality, moulded rubber ferrules
- · Reinforced with metal insert for extended wear
- Product code: BE3652
- UOM: EA

Ferrule, Rubber, Black Type 3/4"

- Latex-free
- Rubber tip flexes for cushioning and shock absorption
- Large surface are provides superior tracking
- High quality, moulded rubber ferrules
- · Reinforced with metal insert for extended wear
- Product code: BE3653
- UOM: EA

Spare Parts / Accessories



Walking Frame Wheels

· A pair of wheels for walking frames

Fits all Breezy Everyday walking frame except BE3637 and BE3638

• Wheel Diameter: 125 mm (5")

• Product Code: BE3624

UOM: Pair

ARTG No: 160745

Breezy Everyday Extension Legs 80mm

• 2 legs with suction cups

• 2 legs with Ferrules

Fits BE98309, BE98310B, BE3645

Product Code: BE3659ARTG No: 160749

Breezy Everyday Extension Legs 80mm

• 2 legs with suction cups

• Fits BE98309, BE98310B, BE3645

Product Code: BE3660

ARTG No: 160749

Breezy Everyday Extension Legs 80mm

2 legs with Ferrules

Fits BE98309, BE98310B, BE3645

Product Code: BE3661ARTG No: 160749

Breezy Everyday Pedal Exerciser, Tension Adjustable Knob

Fits BE26001F

Product Code: BE3662

Breezy Everyday Walker Skis

• Fits 28.6mm shaft

Only fits Ball Walkers

• Product Code: BE3664

ARTG No: 203275

· Product Weight: 24g each

The Breezy Everyday range is imported by:





Sunrise Medical Australia Pty Ltd.

11 Daniel Street Wetherill Park NSW 2164 Australia

Phone: (02) 96786600

Email: orders@sunrisemedical.com.au

www.sunrisemedical.com.au