

**NEW  
PRODUCT**

**GO  
OFFROAD!**

LARGE PNEUMATIC TYRES



COMPACT



# ALL TERRAIN COMPACT SEAT WALKER

RG4412

## FEATURES

- Strong lightweight construction
- User friendly loop brake levers
- Height adjustable handles
- Ergonomic handgrips to reduce pressure
- Brake + park brakes
- Simple brake cable adjustment
- Mesh seating
- Padded Backrest
- Large diameter wheels for excellent mobility
- Pneumatic tyres
- Swivelling front wheels
- Removable bag
- Walking stick holder
- Folds up for transit or storage

## SPECIFICATIONS Dimensions shown are nominal

Overall:

**890mm (L) x 610mm (W) x 900mm (H)**

Overall Folded:

**450mm (L) x 220mm (W) x 670mm (H)**

Handle Height:

**900 - 1025mm (5 settings)**

Seat:

**460mm (W) x 250mm (D) x 600mm (H)**

Wheel Diameter:

**Ø305mm / 12"**

Bag:

**370mm (L) x 200mm (W) x 200mm (D)**

Product Weight:

**11kg**

Maximum User Weight:

**136kg**

Maximum Bag Load:

**2kg**

### IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.



**Mobility**  
Hire & Sales