

Cari Chair [®] carilex

Recliner and Chair Long Overlay

Giving you, the care you need whenever, wherever!



Anti Skid Security



Washable, Waterproof
Cushion Cover



Quick Coupling

Global Leader in Wound Care, Distributed in Over 40 Countries Worldwide

Providing uninterrupted dynamic/static therapy CariChair fits seamlessly with recliner chairs while offering pressure relief and redistribution for patients with sensitive tissues. CariChair Long Overlay is the perfect accessory for patients who need to change from bed to chair or vice-versa.

GENERAL FEATURES

Extra Support at the Sacrum Area	Specially designed deep air cells for excellent immersion and pressure relief
Whole body therapy	Head-to-toe protection for patients even in a sitting position
Seamless Transition from Bed to Chair	Continuous, uninterrupted care whenever, wherever
Quick Release Coupling	Provides an easy connection between power unit and cushion

TECHNICAL DETAILS

Risk/Category	Med - High
Waterlow Score	15+
Pressure Sore Grading	3
Safe Working Load	Up to 150kgs
Air Cell Configuration	26 cells
Cycle Times	10, 15 or 20 minutes
Available Size/Dimensions	180cm x 50cm x 14cm
Travel/Weight	Weighing 2.2kg
Clinical Settings	Suitable for Hospital and Home Care environments
Compatible Pumps	CariChair
Warranty	1 year

carilex[®]

Manufacturing for over 30 Years - 10 Years in the Australian Market