



SELF PROPEL FOLDING MOBILE SHOWER COMMUNE WITH SWINGAWAY FOOTREST - 445MM

Item No. KA130S

Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind . 3 strap back rest for comfort and easy access . Integrated push-handle for manoeuvrability . Swing-up arms for unrestricted side transfer . Swingaway and removable height adjustable footrest . 2x

DESCRIPTION

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind .
- 3 strap back rest for comfort and easy access .
- Integrated push-handle for manoeuvrability .
- Swing-up arms for unrestricted side transfer .
- Swingaway and removable height adjustable footrest .
- 2x 125mm total locking castors .
- 2x 22"self-propel wheels inc. push to lock brakes .
- Designed to be used with a K Care patented removable seats .
- Push rims available as optional extra .
- Padded seat not included .

SPECIFICATIONS

Overall Width	630 mm
Overall Depth	1050 mm
Overall Height	1005 mm
Seat Width	445 mm
Seat Depth	430 mm
Floor To Seat Height	550 mm
Armrest Type	Loop
Backrest Height	400 mm
Footrest Type	Swing Away
Product Weight	15 kg
Max User Weight	200 kg
Brand	K Care