

SB1 ShowerBuddy

For transfers to standard showers

The ShowerBuddy SB1 shower transfer system will turn your existing bathroom into an accessible bathroom that allows you to enjoy showering again without the excessive cost, delays and stress of a full bathroom remodel. With only one transfer, you can once again enjoy a shower in your own bathroom.

The ShowerBuddy SB1 alleviates heavy lifting for family members and caregivers by eliminating multiple transfers during toileting and showering. This unique internationally patented shower transfer system allows the user to toilet and shower in one easy trip.



Product features:

- Weight capacity – 160kg (350lb)
- Rolls over most standard toilets
- Low maintenance
- Arms - removable / lockable / flip-back
- Height adjustable
- Safety mechanisms
- Corrosion resistant Aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Rotates 360 degrees inside shower (subject to space)
- Footrest - fold-up
- Safety belts - lap & chest belts
- Caster Wheels – 125mm (5") x 4 locking
- Commode system
- Tool-free assembly
- Full details see page 31-36


Product Dimensions:

- Rolling Base 467mm W x 560mm D (18.4" x 22")
- Rolling Base including Footrest 467mm W x 675mm D (18.4" x 26.5")
- Width – over arms 570mm (22.5")
- Width – Inside Arms 475mm (18.7") - smaller users see accessory pages
- Width - Inside Main legs 380mm (14.9")
- Back Width 460mm (18")
- Back Rest Height from Cushion 440mm (17")
- Toilet Height Clearance (maximum) 590mm (23.2")
- Toilet clearance can depend upon the shape of the toilet pan
- Bridge Length (Standard) 405mm (16") (or Customized)
- Seat Cushion Dimensions 450 W x 420mm D (17.7" x 16.5")
- Seat Cushion Height from Floor 550mm to 680mm (Adjustable) (21.5" – 26.7" – Adjustable)
- Seat Cushion Opening 175mm W x 335mm L (6.8" x 13")
- Shower Base (Over Suction Cup Feet) 490mm x 490mm (19.2" x 19.2")



SB1 transfers into a standard shower that is not accessible

showerbuddy™



ShowerBuddy™

QUICK REFERENCE GUIDE

FEATURES	Bath Series			Shower							Hybrid	Paediatric	
	SB1	SB2	SB2T	SB3T	SB6c	SB6c-22	SB6c-26	SB6w	SB7c	H5		P1	P2
Height Adjustable	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
Lowest height (mm)	550	560	630	510	500	500	500	570	510	500	Note 1	Note 1	
Lowest height (inches)	21.6"	21.9"	24.8"	20"	19.6"	19.6"	19.6"	22.4"	20"	19.6"	Note 1	Note 1	
Highest height (mm)	620	645	715	600	575	575	575	595	510	575	Note 1	Note 1	
Highest height (inches)	24.4"	25.4"	28"	23.6"	22.6"	22.6"	22.6"	23.4"	20"	22.6"	Note 1	Note 1	
Commode & Toilet Chair	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Tilt N Space	No	No	Yes	Yes	No	No	No	No	No	No	Note 2	Note 2	
Double Locking Struts on Tilt	N.A.	N.A.	Yes	Yes	N.A.	N.A.	N.A.	N.A.	N.A.	N.A.	Note 2	Note 2	
Transfers into Raised Shower	Yes	~	~	No	No	No	No	No	No	Note 7	Yes	Yes	
Transfers into level Shower	~	~	~	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Note 3	Note 3	
Rotates in Shower	Yes	No	No	No	Note 4	Note 4	Note 4	Yes	Note 4	Note 4	Note 5	Note 5	
Transfers over Bath	No	Yes	Yes	No	No	No	No	No	No	Note 7	Note 6	Note 6	
Locking Castor wheels x 4	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	2 Locking	No	N.A.	N.A.	
Neck Rest	Accessory	Accessory	Yes	Yes	No	No	No	No	No	No	Yes	Yes	
Double safety Brakes	Yes	Yes	Yes	N.A.	N.A.	N.A.	N.A.	N.A.	N.A.	Yes	N.A.	N.A.	
Lap / Chest Belts	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	No	Yes	Yes	
Soft Seat Cushion	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Soft Overlay / Day Chair Cushion	Accessory	Accessory	Accessory	Accessory	Accessory	No	No	Accessory	No	No	No	No	
Foldaway Arms	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Removable Arms	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	
Foldaway Footrest	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Removable Footrest	No	No	Yes	Yes	Yes	Yes	Yes	Yes	No	No	Yes	Yes	
Tool less easy assembly	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Non Corrosion Construction	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Weight Capacity Kg	160kg	160kg	160kg	160kg	136kg	170kg	192kg	136kg	120kg	160kg	30kg	90kg	
Weight Capacity Lb	350lb	350lb	350lb	350lb	300lb	375lb	425lb	300lb	265lb	350lb	66lb	200lb	
2 Year Limited Warranty	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	

~

denotes that the chair could be used but would not normally be applied to the situation.

Note 1 P1 and P2 accessory is fitted to SB1, SB2, SB2T or SB3T

Note 2 so the height is dependant upon which model it is fitted to with the seat height matching the model spec + 25mm (1")

Note 3 P1 and P2 can be tilting or non-tilting depending upon which model chair the accessory is fitted to. It tilts when fitted to SB2T and SB3T only.

Note 4 P1 and P2 can be used to transfer into a level shower when using SB3T chair.

Note 5 Could rotate by turning whole chair provided the shower size allowed this.

Note 6 P1 and P2 rotate within the shower when fitted to our SB1 chair only.

Note 7 P1 and P2 will transfer over a bath tub when used in conjunction with SB2 or SB2T chairs

Note 8 H5 Hybrid is designed for those with some mobility, yet requiring medium support. It is anticipated that users could fit the H5 over a shower step or bath tub to use the H5 within the shower or bath tub for support while standing or sitting.