



<11.5 km
Range per Charge



147 kg
Weight Capacity



122.5 cm
Turning Radius



GO GO®
ENDURANCE Li

For people who are always on the GO

The Go Go® Endurance Li is a highly convenient travel companion that provides an effortless mode of transportation for people who are always on the go. This 4-wheel scooter can easily be disassembled into 5 portable pieces, making it easy to transport or store in the trunk of a car. With the lightweight Lithium ion battery pack, you can travel up to 23km on a single charge, giving you plenty of freedom to explore your surroundings and reach more destinations.



GO GO[®] ENDURANCE Li



Long-lasting, lithium-ion battery



Feather-touch disassembly




CTS Suspension on front and rear

FRONT TYRES ²	5 x 20.3 cm
REAR TYRES ²	6.4 x 18 cm
MAXIMUM SPEED ¹	6.7 km/h
GROUND CLEARANCE ²	3.6 cm at motor
TURNING RADIUS ²	122.5 cm
LENGTH & WIDTH ²	L: 107 cm W: 52 cm
SEAT-TO-GROUND HEIGHT ²	53.9 - 57.7 cm
STANDARD SEAT	Compact, Foldable L: 43.1 cm (17") W: 43.1 cm (17") Weight: 10.4 kg L: 50.8 cm (20") W: 45.7 cm (18") Weight: 14.1 kg
TOTAL WEIGHT WITHOUT BATTERIES ²	46.5 kg

SEAT-TO-DECK HEIGHT ²	41.2 cm - 45 cm
BATTERY WEIGHT	3.6 kg
BATTERY REQUIREMENTS ²	(2) 8AH
MAXIMUM RANGE PER CHARGE ^{1,3}	Up to 11.5 km
BATTERY CHARGER	Off-board, 3.5 amp (Lithium Ion)
SUSPENSION	CTS Front and Rear
MAXIMUM WEIGHT CAPACITY	147 kg
STANDARD BODY COLOURS	■ Granet Red ■ Euro Grey
WEIGHT OF HEAVIEST PIECE ²	19 kg (front section)
3W/4W OPTION	4W Only

1. Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of (+/-) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of (+/-) 3%.
3. Results derived from theoretical calculations based on lithium battery specifications and drive system performance. Test conducted at maximum weight capacity.

 Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go Go[®] Endurance Li and its lithium ion battery are permitted.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.



PRIDE[®]
MOBILITY
LIVE YOUR BEST[™]

20-24 Apollo Drive, Hallam, Victoria 3803 | (03) 8770 9600 | www.pridemobility.com.au

