

The Invacare Colibri is a stylish, simple and colourful micro lightweight scooter, designed for those who enjoy a totally independent lifestyle. Less than a minute is all that's needed to disassemble the Colibri for transportation.





# Colibri

Thanks to the unique patented Invacare LiteLock™ system the scooter can easily be taken apart without the need for tools and fits neatly into any car boot. Plus, with its small footprint the Invacare Colibri offers greater manoeuvrability in tight, or limited spaces.

Additionally, you can choose from a range of trendy colour panels to match your style, swapping them according to your needs.

Choose where to charge your battery: directly on the Colibri or detached from it.

The soft, padded seat provides comfort during driving. For increased comfort there is also the optional 'easy fit' seat suspension kit.

#### **Features**



Effortless dismantling can be done in less than a minute for transportation.



A soft padded seat provides additional comfort

### **Tech Specs**



Backrest width 465 mm



Seat depth 400 mm



Floor to seat height 370 - 420 mm



**Backrest height** 385 mm



Driving unit width



Total length 1010 mm



**Empty weight** 49 kg (18 Ah)



Max user weight 136 kg



**Battery capacity** 2 x 18 Ah



Turning diameter 2200 mm



Max climable obstacle height 45 mm



Ground clearance 50mm



Max safe slope 6° / 10.5°



Max driving distance Up to 16 km



Max driving speed 8 km/h



Motor capacity 200 W

## **Colour Options**



Jade Green

Sapphire Blue



Diamond Silver



Marble White



#### FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au Asia: Phone +61 2 8839 5330, email anzexportsales@invacare.com or visit www.invacare.com Thailand: Phone +66 0 2821 5515 or visit www.invacare.com



