Pegasus Pro

COMFORT, SAFETY AND RELIABILITY WITH ADDED PERFORMANCE



With enhanced battery performance, superior suspension and large 12" wheels, the Pegasus PRO is ideal for users who need that extra performance. comfort, safety and reliability features to make it the ideal choice for users who demand more.



Making Life's Experiences Possible



Pegasus Pro

The advanced suspension system on the Pegasus PRO has been designed to bring better performance and give a smoother ride over a variety of surfaces and terrains. In addition, the Pegasus PRO features a powerful motor and 12" pneumatic tyres to ensure that whether you are in the city or the country, you get where you need to go in comfort and style.

An on-board electronic speed reduction system automatically reduces speed around bends, ensuring users are confidently making safe, comfortable manoeuvres. Furthermore, a two-step disengaging lever prevents the scooter from freewheeling if it is knocked out of drive mode.

The electronics and motor have additional protection to guard against water and corrosion, safeguarding the vital workings of the scooter against any unnecessary malfunctions.

Features



Users can easily adjust the tiller with a lever to suit their needs.



Clear, bright brake lighting warns people that the scooter is slowing down.



Ergonomically designed steering system makes for easier, more responsive control.



Dashboard configuration.



Ensures immediate braking if required.

Tech Specs



510 mm

Battery capacity

2 x 12 V/75 Ah



Seat depth 470 mm

Turning diameter

2750 mm



Max climable

obstacle height

100 mm

Backrest height 490mm



100 mm



10



Max driving distance 52 km

Total length

1320 mm



136 kg

Max driving speed

10 km/h

KM/H



Refer to order form for product codes

Max user weight 160 kg



Motor capacity 550 W / 1300 W

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au Asia: Phone +61 2 8839 5330, email anzexportsales@invacare.com or visit www.invacare.com Thailand: Phone +66 0 2821 5515 or visit www.invacare.com



