



COMFORTLITE TRANSIT WHEELCHAIR

RG312GR18: GREEN RG312BL18: STEEL BLUE

FEATURES

Lightweight: 9.5kg (Aluminium frame) Fixed armrest with PU padding Foldable backrest with storage pocket Removable seat and backrest cushions Double Handbrake

Fixed leg rests with flip up footrests Large flat free tyres

Seat belt

Calf strap (removable) / Side straps (removable) Anti tippers

SPECIFICATIONS (Dimensions shown are nominal)

Overall:

840mm (L) - Footrests up 610mm (W) - 270mm (Folded) 975mm (H) - 730mm (Folded)

480mm (H) x 460mm (W) x 420mm (D)

Backrest:

450mm (H) x 445mm (W)

Armrest:

200mm (H)

Wheels:

150mm / 6" (Front) 300mm / 12" (Rear)

Product Weight:

9.5kg

Max. User Weight:

100kg

Please contact **REDGUM** if you require more information on this product.



P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090 Part of the FOR-DE Group Pty Ltd.

redgumbrand.com.au

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.